



## Cylinder, graduated, double metric, 1 L

**Proizvod #:** 2928305  
**EUR Cijena:** Kontaktirajte nas  
**Nije prikazan datum isporuke**

These glass pyrex@ graduated cylinders feature two sets of graduations for reading either up or down Calibrated "to deliver"

---

### Specifikacije

|                           |                    |
|---------------------------|--------------------|
| Certificiranja:           | Class A            |
| Kapacitet:                | 1000 mL            |
| Materijal tijela mjerača: | Glass              |
| Mjerna metoda:            | To Deliver         |
| Visina:                   | 465 mm             |
| Volumen/ veličina paketa: | 1 each             |
| Vrsta:                    | Dvostruko metrički |