



## Bottle, BOD, 300 mL, 6/pk

|                              |                   |
|------------------------------|-------------------|
| <b>Proizvod #:</b>           | <b>62106</b>      |
| EUR Cijena:                  | Kontaktirajte nas |
| Nije prikazan datum isporuke |                   |

Set of 6 glass BOD bottles, 300 mL, with stoppers. For use with LBOD sensor to determine the Biological Oxygen Demand. Borosilicate clear glass bottle with flared lip and ground-glass stopper to prevent air entrapment and for water seal. Frosted marking area. Stoppers included.

---

### Specifikacije

|                           |                    |
|---------------------------|--------------------|
| Body Material:            | Glass              |
| Broj metode:              | 8043               |
| Kapacitet:                | 300 mL             |
| Materijal:                | Borosilicate Glass |
| Package size UOM:         | 6pk                |
| Volumen/ veličina paketa: | 6 /pk              |