



Aplikacijski softver za enzimsku analizu hrane

Proizvod #:

LZW943

EUR Cijena:

Kontaktirajte nas

Dostupno

The application software LZW937 is a compilation of enzymology tests

manufactured by R-Biopharm AG, Darmstadt that can be performed with the spectrophotometer DR6000 and carousel insert LZW902.99.00001. The methods can be modified for individual settings.

Methods are ready to go.

The methods are convenient to upload via USB.

Keep control of time depending multiple measurements.

A timer function alerts about the next step and starts evaluation of multiple samples automatically.

Get the final result immediately after the analysis.

The photometer does the complex result calculation for you.

Potrebni dodaci

- Rotirajući držač 7 x 1 cm (Item LZW902.99.00002)