



Be Right™



POCKET Colorimeter II Colorimeter Test Kit for Manganese (HR) analysis

Proizvod #: 5870015

Opasno po zdravlje

NEDOSTUPAN ARTIKL
Ovaj artikl više nije dostupan.



Small in size, big on waterproof analysis

Portable Colorimeter programmed for the determination of Manganese, 0.20-20.0 mg/L Mn, Periodat Oxidation Method. In carrying case, complete with reagents, cuvettes and operating instructions.

Simple: all functions are available via only four keys

Power: battery operation for a maximum of 2000 tests

Clear readings: even under difficult conditions, thanks to backlit display and large readout

Reliable results with no mains connection

Rugged construction yet lightweight

Specifikacije

Broj metode:	8034
Broj testova:	100
Compliance certifications:	CE
Detektor:	Silicon detector
Dimenzije (V x Š x D):	155 mm x 61 mm x 35 mm
Garancija:	2 godine
Klasa kućišta:	IP67 (isključuje odjeljak za baterije)
Masa:	0,23 kg
Mjerni rang:	0,2 - 20,0 mg/L Mn
Naziv metode:	Cold periodate
Parametar:	Manganese
Smallest increment:	0.1
Valna duljina:	Fixed wavelength \pm 2 nm varies with model
Vrsta kućišta:	D
Zaslon:	LCD, backlit