



Standard za nitrat, 400 mg/L NO₃ (90,4 mg/L NO₃-N), 500 mL

Proizvod #:

LCW863

EUR Cijena:

Kontaktirajte nas

Dostupno

Standardna otopina za online određivanje nitrata

Specifikacije

Koncentracija:	400 mg/L NO ₃ (90,4 mg/L NO ₃ -N)
Parametar:	Nitrat
Platforma:	Otopina
Uređaj:	Nitratax
Volumen/ veličina paketa:	500 mL