



Reagent powder pillows, FerroVer iron, 10 mL, pk/1000

Proizvod #:

2105728

EUR Cijena:

Kontaktirajte nas

Dostupno

Opasno po zdravlje



Powder reagent in sealed package, pk/1000. For 10mL sample.

Specifikacije

Broj metode:	8008
Broj testova:	1000
Čuvanje uređaja:	10 °C - 25 °C
EPA compliant:	Da
Mjerni rang:	0,02 - 3,00 mg/L Fe
Naziv metode:	FerroVer
Parametar:	Iron
Platforma:	Praškasti reagensi PP
Uređaj:	All except DR/3000, DR/2000
Volument/ veličina paketa:	1000 /pk