



Molybdenum 1 Reagent Powder Pillows, 20 mL, pk/100

Proizvod #:	2352449
EUR Cijena:	Kontaktirajte nas
Dostupno	

Powder reagent in sealed package for the determination of Molybdenum by the ternary complex method (Hach method 8169), 0.02 - 3.00 mg/L Mo. For 20 mL sample.

Specifikacije

Broj testova:	100
Naziv metode:	Ternary Complex
Platforma:	Praškasti reagensi PP
Volumen/ veličina paketa:	100 /pk