



CuVer 1 Copper Reagent Powder Pillows, 0.04-5.00mg/L Cu

Proizvod #:

1418899

Opasno po zdravlje

EUR Cijena:

Kontaktirajte nas

Dostupno



Powder reagent in sealed package for the determination of Copper, Bicinchoninate Method. For 25mL sample. 100/pack

Specifikacije

Broj metode:	8506
Broj testova:	100
EPA compliant:	Da
Mjerni rang:	0,04 - 5,00 mg/L Cu
Naziv metode:	Bicinchoninate
Parametar:	Copper
Platforma:	Praškasti reagensi PP
Uređaj:	DR/3000, DR/2000
Volumen/ veličina paketa:	100 /pk