



PAN Indicator Solution, 0.1% (50mL SCDB)

Proizvod #:

2122426

EUR Cijena:

Kontaktirajte nas

Dostupno

Opasno po zdravlje



Colorimetric reagent for low level Manganese determinations. Self-contained dropping bottle.

Specifikacije

Broj metode:	8149
Mjerni rang:	0,006 - 0,700 mg/L Mn
Naziv metode:	PAN
Parametar:	Manganese
Platforma:	Otopina
Volument/ veličina paketa:	50 mL