



Porphyrin 2 Reagent Powder Pillows, 25 mL, pk/100

Proizvod #:

2187569

Opasno po zdravlje

EUR Cijena:

Kontaktirajte nas

Otprema unutar 2 tjedna



Powder reagent in sealed package for the determination of Copper by the Porphyrin method, 2 - 210 µg/L Cu. For 25 mL sample.

Specifikacije

Naziv metode:

Porphyrin

Platforma:

Praškasti reagensi PP

Volumen/ veličina paketa:

100 /pk