



Cobalt / Nickel Reagent Set

Proizvod #:

2242600

EUR Cijena:

Kontaktirajte nas

Dostupno

Opasno po zdravlje



For the determination of Cobalt, 0.01-2.00 mg/L Co, and Nickel, 0.006-1.000 mg/L Ni, PAN Method. For 25 mL sample.

Specifikacije

Broj metode:	8150
Broj testova:	100
Mjerni rang:	0,01 - 2,00 mg/L Co
Mjerni rang (2):	0,006 - 1,000 mg/L Ni
Naziv metode:	PAN
Opis:	Nickel Reagent Set, PAN Method, 25 mL
Parametar:	Cobalt, Nickel
Posebna značajka:	Cobalt Nickel Reagent set, 25ml
Uređaj:	DR/3000, DR/2000, DR900
Volumen uzorka:	25 mL