



Reagent set, copper, porphyrin, 25 mL

Proizvod #:

2242700

Opasno po zdravlje

EUR Cijena:

Kontaktirajte nas

Nije prikazan datum isporuke



Powder reagents in sealed packages for the determination of Copper, Porphyrin Method. For 25 mL sample.

For determination of copper by the Porphyrin Method. Hach Method 8143. Range: 2 - 210 µg/L Sample Size: 25 mL. Reagent Set includes: Copper Masking Reagent Powder Pillows, Porphyrin 1 and 2 Reagent Powder Pillows for 100 tests.

Specifikacije

Broj metode:	8143
Broj testova:	100
Mjerni rang:	2 - 210 µg/L
Naziv metode:	Porphyrin
Platforma:	Praškasti reagensi PP
Uređaj:	DR/3000, DR/2000