



QAC reagent 1, 0.2- 5.0 mg/L CTAB, 50 powder pillows

Proizvod #: 2401066
EUR Cijena: Kontaktirajte nas
Nije prikazan datum isporuke

Powder reagent in sealed package for the determination of Quaternary Ammonium Compounds (QAC).

Specifikacije

Broj metode:	8336
Broj testova:	50
Naziv metode:	Direct Binary Complex
Parametar:	Quaternary Ammonium Compounds
Platforma:	Praškasti reagensi PP
Volumen/ veličina paketa:	50/pk