



Klor, tekući DPD komplet reagensa, 0,02-2,00 mg/L Cl₂

Proizvod #:

HPT310

Opasno po zdravlje

EUR Cijena:

Kontaktirajte nas

Dostupno



Specifikacije

Broj testova:	100
Čuvanje uređaja:	15 °C - 25 °C
Mjerni rang:	0,02 - 2,00 mg/L Cl ₂
Naziv metode:	DPD
Opis:	Slobodni/ukupni klor
Parametar:	Klor
Platforma:	Otopina
Posebne napomene:	za procjenu su potrebne kružne bočice od 1 inča
Standard method:	ISO 7393-2