



Free Chlorine Reagent Powder Pillows, DPD (pk/100)

Proizvod #:

1407799

Opasno po zdravlje

EUR Cijena:

Kontaktirajte nas

Dostupno



For determination of free chlorine by the DPD method. For use with Hach's Chlorine Visual Test Kits.

Sample Size: 5 mL. Range: up to 3.5 mg/L Cl_2 .

Pack of 100 powder pillows.

Boric Acid free!

Specifikacije

Broj testova:	100
Čuvanje uređaja:	10 - 25 °C
Mjerni rang:	5 mL
Naziv metode:	DPD
Parametar:	Chlorine, Free
Platforma:	Praškasti reagensi PP