



Total Hardness Reagent Set for Digital Titrator (LR)

Proizvod #:

2448000

EUR Cijena:

Kontaktirajte nas

Otprema unutar 2 tjedna

Opasno po zdravlje



For the determination of Total Hardness, 10-160 mg/L CaCO₃, EDTA Method.

Specifikacije

Broj metode:	8213
Broj testova:	40-100
EPA compliant:	Da
Mjerni rang:	10 - 160 mg/L CaCO ₃
Naziv metode:	EDTA Titration
Parametar:	Hardness (Total)
Platforma:	Digitalni Titrator
Posebna značajka:	Reagent set, LR, Digital Titrator