



EDTA solution, 1M, 50mL SCDB

Proizvod #:

2241926

Opasno po zdravlje

EUR Cijena:

Kontaktirajte nas

Dostupno



For determination of calcium and magnesium hardness in water by the Calmagite Colorimetric method, Hach Method 8030. Also used for the determination of calcium and magnesium in ultrapure saturated NaCl feed brine (HNB method). 50 mL self-contained Dropping Bottle.

Specifikacije

Broj metode:	8030
Naziv metode:	Calmagite Colorimetric
Parametar:	Calcium and Magnesium Hardness
Platforma:	Otopina
Volumen/ veličina paketa:	50 mL

Potrebni dodaci

- Zaštitne rukavice veličine L, plave, nitril, bez pudera, 50 komada (Item SM743L)