



Komplet reagensa za fenole (bez dostave avionom)

Proizvod #:

2243900

Opasno po zdravlje

EUR Cijena:

Kontaktirajte nas

Dostupno



Za određivanje fenola, 4-aminoantipirin metoda.

Specifikacije

Broj metode:	8047
Broj testova:	100
EPA compliant:	Da
Mjerni rang:	0,002 - 0,200 mg/L Fenol
Naziv metode:	4-aminoantipirin
Opis:	Komplet reagensa za fenole
Parametar:	Fenoli
Shipping Method:	No Air Shipments
Uređaj:	Sve osim DR/800, PC II