



## Standardna otopina fosfata, 15 mg/L kao $\text{PO}_4$ (NIST), 100 mL

**Proizvod #:** 1424342  
**EUR Cijena:** Kontaktirajte nas  
**Otprema unutar 2 tjedna**

For accuracy checks in phosphate (phosphorus) determinations

---

### Specifikacije

|                           |                                                           |
|---------------------------|-----------------------------------------------------------|
| Koncentracija:            | 15 mg/L $\text{PO}_4$                                     |
| Opis:                     | Standardna otopina fosfata, 15 mg/L $\text{PO}_4$ , 100mL |
| Parametar:                | Fosfat                                                    |
| Platforma:                | Otopina                                                   |
| Volumen/ veličina paketa: | 100 mL 100                                                |