



Be Right™



1990 $\mu\text{S}/\text{cm}$ standard za vodljivost; (995 mg/L TDS, 1.000 mg/L NaCl), 1L

Proizvod #: 210553
HRK Cijena: Kontaktirajte nas
Poziv za datum slanja

1.000 mg/L NaCl. Za sonde za vodljivost.

Prikladno i za standardiziranje sklopa za vodljivost, $\mu\text{S}/\text{cm}$ ($\mu\text{mhos}/\text{cm}$).

Specifikacije

Koncentracija:	1990 $\mu\text{S}/\text{cm}$
Opis:	Standardna otopina natrijevog klorida, 1.990 $\mu\text{S}/\text{cm}$ (NIST), (995 mg/L TDS, 1.000 mg/L NaCl), 1 L
Parametar:	Vodljivost
Platforma:	Otopina
Volument/ veličina paketa:	1 L