



Standard Solution, chlorine, 50-75 mg/L as Cl₂, pk/20 - 2 mL PourRite ampules

Proizvod #:	1426820
EUR Cijena:	Kontaktirajte nas
Dostupno	

For accuracy checks in chlorine determinations. Pack of 20 2-mL ampules. Lot-specific concentration displayed on label.

Simple Accuracy Checks

Voluette and PourRite Ampule Standards are designed to work directly with Standard Additions accuracy checks in Hach procedures and photometer software. Use the TenSette Pipet to add precise spikes of standard and check percent recovery, to confirm analytical system performance.

Minimize waste

Open a fresh standard ampule for each accuracy check. Each ampule contains enough standard for one set of standard additions.

Reliable results

Sealed ampules prevent accidental contamination or evaporation of standard solutions.

Specifikacije

Koncentracija:	50 mg/L – 75 mg/L Cl ₂
Platforma:	Otopina