



Otopina nitrita, zaliha, 250 µg/mL (mg/L) kao N, 500 mL

Proizvod #:

2340249

Opasno po zdravlje

EUR Cijena:

Kontaktirajte nas

Otprema unutar 2 tjedna



For accuracy checks in Nitrite determinations.

Za provjeru preciznosti u određivanju nitrita. Boca od 500 mL.

Specifikacije

Koncentracija:	250 mg/L N
Opis:	Zaliha otopine nitrita, 250 µg/mL kao N, 500 mL, APHA
Parametar:	Nitrit
Platforma:	Otopina
Volumen/ veličina paketa:	500 mL 500