



Stablcal Stabilised Formazine Turbidity Standard, 20 NTU, 100 mL

Proizvod #:

2660142

Opasno po zdravlje

EUR Cijena:

Kontaktirajte nas

Dostupno



Ready for use, stabilised true Formazine Primary standard solution, 20 NTU, 100 mL bottle. Suited for all turbidimeters. USEPA accepted for calibration purposes. For TL23, TU52, 2100Q, 2100N, 2100AN series.

Precise, Pre-mixed Formazin Solutions

With proprietary manufacturing technology, Hach prepares Stablcal Standards in precise concentrations at a wide variety of concentrations - providing accurate turbidimeter calibrations.

Based on Formazin Primary Standards

Stablcal Stabilized Formazin Standards contain the same light scattering polymer as traditional Formazin primary turbidity standards - and perform equivalently.

Enhanced Stability Saves Time and Labor

While Formazin primary standards can be prepared directly in the lab, the process is labor intensive and time consuming - and the diluted standards are not stable. Stablcal is a stabilized version of Formazin - providing repeatable results, without the time and effort required for lab-prepared standards.

Confidence in results

Hach was the first to apply Formazin as a turbidity standard in the drinking water industry, more than 40 years ago. Stablcal standards are based on decades of experience in turbidity measurement.

Specifikacije

Koncentracija: 20 NTU

Platforma: Otopina

Potrebni dodaci

- Zaštitne rukavice veličine L, plave, nitril, bez pudera, 50 komada (Item SM743L)