



NitriVer 3 Nitrite Reagent Powder Pillows, 0.02-0.30 mg/L NO₂⁻-N

Proizvod #:

1406599

Opasno po zdravlje

EUR Cijena:

Kontaktirajte nas

Dostupno



Powder reagent in sealed package for the determination of Nitrogen, Nitrite. 25mL sample, pk/100

Specifikacije

Broj metode:	8507
Broj testova:	100
Mjerni rang:	0,002 - 0,300 mg/L NO ₂ ⁻ -N
Naziv metode:	Diazotization
Platforma:	Praškasti reagensi PP
Uređaj:	DR/2000, DR/3000
Volumen/ veličina paketa:	100 /pk