



Be Right™



NitriVer 3 Nitrite Reagent Powder Pillows, 5 mL, pk/100

Proizvod #:

1407899

Opasno po zdravlje

EUR Cijena:

Kontaktirajte nas

Dostupno



For nitrite determination by the NitriVer 3 diazotization method. Sample size 5 mL. For use with Hach visual test kits. Range: up to 0.5 mg/L NO₂-N.

Specifikacije

Naziv metode:

Diazotization

Platforma:

Praškasti reagensi PP

Potrebni dodaci

- Zaštitne rukavice veličine L, plave, nitril, bez pudera, 50 komada (Item SM743L)