



SPADNS 2 Fluoride Reagent, arsenic-free, ACCUVAC, 0.02-2.00 mg/L F

Proizvod #:

2527025

EUR Cijena:

Kontaktirajte nas

Dostupno

Opasno po zdravlje



ACCUVAC ampoules contain the precise amount of reagent for a single test and can be used as a measurement cuvette. For the determination of Fluoride by the SPADNS 2 method. pk/25

Specifikacije

Broj testova:	25
EPA compliant:	Da
Method Number 1 Chemistry:	10225
Mjerni rang:	0,02 - 2,00 mg/L F ⁻
Naziv metode:	SPADNS 2
Parametar:	Fluoride
Platforma:	Accuvac
Volumen/ veličina paketa:	25 /pk