



## SPADNS 2 (Arsenic-free) Fluoride Reagent Solution, 500 mL

**Proizvod #:**

**2947549**

EUR Cijena:

Kontaktirajte nas

Dostupno

Opasno po zdravlje



For determination of Fluoride by Hach SPADNS 2 Method 10225. The SPADNS2 reagent formulation removes interference from chlorine without arsenic. Range: 0.02 - 2.00 mg/L F<sup>-</sup>. EPA compliant. 500 mL bottle

### Identical performance to traditional Hach SPADNS

Use the same Hach instruments, test procedures, and calibration curves

---

## Specifikacije

|                |                                 |
|----------------|---------------------------------|
| Broj metode:   | 10225                           |
| Broj testova:  | 125                             |
| EPA compliant: | Da                              |
| Mjerni rang:   | 0,02 - 2,00 mg/L F <sup>-</sup> |
| Naziv metode:  | SPADNS 2                        |
| Platforma:     | Otopina                         |