



## Fosfor, reaktivni, komplet reagensa, aminokiseline

**Proizvod #:**

**2244100**

**Opasno po zdravlje**

EUR Cijena:

Kontaktirajte nas

Dostupno



Za određivanje fosfata metodom aminokiselina.

---

### Specifikacije

|               |                                   |
|---------------|-----------------------------------|
| Broj metode:  | 8178                              |
| Broj testova: | 100                               |
| Mjerni rang:  | 0,23 - 30,00 mg/L PO <sub>4</sub> |
| Naziv metode: | Aminokiseline                     |
| Parametar:    | Fosfor                            |
| Platforma:    | Otopina                           |
| Uređaj:       | Sve osim PC II                    |