



Reagens za molibdovanadat, 500 mL

Proizvod #:

2076049

EUR Cijena:

Kontaktirajte nas

Dostupno

Opasno po zdravlje



Za određivanje orto fosfata, 0,3-45,0 mg/L PO₄

Specifikacije

Broj metode:	8114
Broj testova:	250
Mjerni rang:	0,3 - 45,0 mg/L PO ₄
Naziv metode:	Molibdovanadat
Parametar:	Fosfor
Platforma:	Otopina
Uređaj:	Sve osim DR/850, DR/820, PC II
Volument/ veličina paketa:	500 mL