



Chlorine Standard, 100-130 mg / L 20 mL Amp

Proizvod #:	2835801
EUR Cijena:	Kontaktirajte nas
Dostupno	

This chlorine standard can be diluted down to verify the accuracy of the CL17 chlorine analyser or create a separate calibration curve. Under normal circumstances, a separate calibration curve is not required for the CL17. Package of one 20 mL ampule.

Specifikacije

Čuvanje uređaja:	Keep refrigerated.
Mjerni rang:	N/A
Uređaj:	CL17

Potrebni dodaci

- Zaštitne rukavice veličine L, plave, nitril, bez pudera, 50 komada (Item SM743L)