



Stablcal Turbidity standard solution 1.0 NTU

Proizvod #:

2659853

Opasno po zdravlje

EUR Cijena:

Kontaktirajte nas

Dostupno



Stabilised Formazin turbidity standard for use in any turbidimeter. 1000mL bottle.

Specifikacije

Čuvanje uređaja:

Store between 5 - 25 °C, away from direct sunlight

Koncentracija:

1.0 NTU

Mjerni rang:

N/A

Parametar:

Turbidity

Volument/ veličina paketa:

1 L