




EDTA Digital Titrator Cartridge, 0.1428 M

Proizvod #:

1496001

EUR Cijena:

Kontaktirajte nas

Nije prikazan datum isporuke

Opasno po zdravlje



Ready-to-use titrant solution for the analysis of Total Hardness, 1-16 G.d.H. CaCO₃.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifikacije

Broj testova:	40 - 100
Mjerni rang:	1 - 16 G,d,h,
Naziv metode:	EDTA (CalVer 2 und ManVer 2)
Parametar:	Hardness, Total
Platforma:	Digitalni Titrator
Titrant:	0.1428 M EDTA

Potrebni dodaci

- Zaštitne rukavice veličine L, plave, nitril, bez pudera, 50 komada (Item SM743L)