



Ruksak za prijenosne mjerače, velik, najlon, kotači, sa pretincima

Proizvod #:

2990800

EUR Cijena:

Kontaktirajte nas

Otprema unutar 2 tjedna

Za udobno prenošenje prijenosnih mjerača, dodataka, i standardnih otopina.

Prikladan

The answer to hands-free portability of analytical testing equipment, these backpacks combine contemporary styling with innovative design.

Prilagodljiv

The adjustable divider system allows you to create customized compartments to organize equipment just the way you like it.

Ergonomski

Lightweight, thick padding material gives superior protection during transportation. A stitched foam back panel makes the backpack well suited for either short or long testing adventures. The sternum buckle system adds support when carrying heavier equipment.

A telescoping handle system and a comfort pad provide additional comfort when used as a backpack. The large backpack can also be used as a wheeled bag or a backpack with in-line skate wheels.

Specifikacije

Dimenzije (V x Š x D):	432 mm x 356 mm x 152 mm
Duljina:	152 mm
Masa:	Kapacitet: 13,6 kg
Materijal:	Najlon, plastični tekstil