



Standardna otopina natrijevog klorida, 100 μ S/cm, 500 mL

Proizvod #: 2971849
EUR Cijena: Kontaktirajte nas
Dostupno

Conductivity standard solution, 100 μ S/cm $\pm$ 1 μ S/cm at 25°C, 49 mg/L $\pm$ 0.25 mg/L as Sodium Chloride (NaCl), (1 μ S/cm $\sim$ 0.5 mg/L TDS), 500 mL bottle.

Specifikacije

Čuvanje uređaja:	10 - 25 °C (protect from light)
Koncentracija:	100 μ S/cm at 25°C $\pm$ 1 μ S/cm at 25°C
Opis:	Natrijev klorid, 100 μ S/cm, 500 mL
Package settings:	Bottle
Parametar:	Vodljivost
Platforma:	Otopina
Posebne napomene:	Observed shelf life after opening, following recommended storage and use: 3 months
Volument/ veličina paketa:	500 mL