



Standardna otopina natrijevog klorida, 10.000 μ S/cm, 500 mL

Proizvod #: 2972249
EUR Cijena: Kontaktirajte nas
Otprema unutar 2 tjedna

Conductivity standard solutions for laboratory and field analyses

Conductivity standard solution, 10000 μ S/cm $\pm$ 100 μ S/cm at 25°C, 5.556 g/L $\pm$ 0.03 g/L as Sodium Chloride (NaCl), (1 μ S/cm $\sim$ 0.5 mg/L TDS), 500 mL bottle.

Traceable to NIST standard reference materials

Whether in the laboratory or in the field, always perform measurement results you can trust is essential. The quality of the standard solutions that are used are key.

Design for use with Hach meters portfolio and other meters.

The basis for accurate conductivity measurements is good care and calibration of the conductivity cell.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Specifikacije

Čuvanje uređaja:	10 - 25 °C (protect from light)
Koncentracija:	10000 μ S/cm at 25°C $\pm$ 100 μ S/cm at 25°C
Opis:	Natrijev klorid, 10.000 μ S/cm, 500 mL
Package settings:	Bottle
Parametar:	Vodljivost
Platforma:	Otopina
Volumen/ veličina paketa:	500 mL