



## Filters, Membrane, 47-mm, 0.45- $\mu$ m, gridded, sterile, Millipore

**Proizvod #:**

**2936100**

EUR Cijena:

Kontaktirajte nas

Nije prikazan datum isporuke

Sterile membrane filters for use in microbiology membrane filtration methods.

---

### Specifikacije

|                           |                       |
|---------------------------|-----------------------|
| Materijal:                | Mixed-Cellulose Ester |
| Promjer:                  | 47 mm mm              |
| Veličina pora:            | 0.45 $\mu$ m          |
| Volumen/ veličina paketa: | 150                   |