



Filter paper, glass fiber, extra thick, pore size 1 µm, diameter 47 mm, 100/pk

Proizvod #:

2842100

EUR Cijena:

Kontaktirajte nas

Nije prikazan datum isporuke

For filtration of heavily contaminated samples. High particulate-holding capacity makes discs efficient depth filters and allows for filtration of large volumes of solutions. Pall Life Sciences.

Specifikacije

Materijal:	Glass Microfiber
Promjer:	47 mm mm
Veličina pora:	1 µm
Volumen/ veličina paketa:	100