



Filter membrane, EZ-Pak, 0.45 micron, 47 mm, black, gridded, 600/pk (millipore)

Proizvod #:	2915900
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Dispenser provides superior convenience and assures sterility. Unique one press on the lever and a sterile membrane is presented (horizontally and inverted) for efficient and convenient assembly. The EZ Pak membranes come as a band of membranes that easily loads into the dispenser. Packaging materials are wound into the dispenser ensuring a clear operating area and quick change of bands. The membranes are Type HA, which signifies compliance with Standard Methods criteria.

Specifikacije

Materijal:	Mixed-Cellulose Ester
Promjer:	47 mm mm
Veličina pora:	0.45 μ m
Volumen/ veličina paketa:	600