



Titration Application Pack Water: pH & Alkalinity

Proizvod #: AP0001.AT1102

EUR Cijena: Kontaktirajte nas

Otprema unutar 2 tjedna

Application Pack for the determination of pH and Alkalinity in water. Suited for Titrilab AT1000 Series, Automatic Titrator Models AT1102/AT1112/AT1122/AT1222.

Specifikacije

Garancija: 24 mjeseci

Kompatibilni instrumenti: AT1102

AT1112

AT1122

AT1222

Parametar: pH/Alkalinity

Primjena: Municipal Applications

Što je u pretincu: Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper.

Vrsta: Application Pack

Što je u pretincu

Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper.