



Titration Application Pack Water: pH, Alkalinity & Conductivity

Proizvod #: AP0002.AT1102
EUR Cijena: Kontaktirajte nas
Otprema unutar 2 tjedna

Application Pack for the determination of pH, Alkalinity and Conductivity in water. Suited for Titralab AT1000 Series, Automatic Titrator Models AT1102/AT1112/AT1122/AT1222.

Specifikacije

Garancija:	24 mjeseci
Kompatibilni instrumenti:	AT1102
	AT1112
	AT1122
	AT1222
Parametar:	pH/Alkalinity & Conductivity
Primjena:	Municipal Applications
Što je u pretincu:	Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper.
Vrsta:	Application Pack

Što je u pretincu

Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper.