



Titration Application Pack Water: Ca & Mg Hardness (ISE)

Proizvod #: AP0003.AT1112

EUR Cijena: Kontaktirajte nas

Nije prikazan datum isporuke

Application Pack for the determination of Calcium and Magnesium Hardness in water. Suited for TitrLab AT1000 Series, Automatic Titrator Models AT1112/AT1122/AT1222.

Specifikacije

Garancija: 24 mjeseci

Kompatibilni instrumenti: AT1112

AT1122

AT1222

Parametar: Ca & Mg Hardness (ISE)

Primjena: Municipal Applications

Što je u pretincu: Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper.

Vrsta: Application Pack

Što je u pretincu

Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper.