



Titration Application Pack, Food & Beverage: pH, Total Acidity & SO₂ in wine

Proizvod #: AP0013.AT1222

EUR Cijena: Kontaktirajte nas

Otprema unutar 2 tjedna

Application Pack for the determination of pH, Total Acidity and Free & Total SO₂ in wine. Suited for Titrator Model AT1000 Series, Automatic Titrator Model AT1222.

Specifikacije

Garancija:	24 mjeseci
Kompatibilni instrumenti:	AT1222
Parametar:	pH, Total Acidity, Free & Total SO ₂ in Wine
Primjena:	Food & Beverage Applications
Što je u pretincu:	Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper, 1 L glass bottle.
Vrsta:	Application Pack

Što je u pretincu

Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper, 1 L glass bottle.