



Titration Application Pack Karl Fischer: Moisture Content

Proizvod #: AP0014.KF1121

EUR Cijena: Kontaktirajte nas

Otprema unutar 2 tjedna

Application Pack for the determination of Moisture content (Karl Fischer). Suited for Titrallab KF1000 Series, Automatic Titrator Model KF1121.

Specifikacije

Garancija:	24 mjeseci
Kompatibilni instrumenti:	KF1121
Parametar:	Moisture Content (Karl Fischer)
Primjena:	Moisture Content Application
Što je u pretincu:	Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper.
Vrsta:	Application Pack

Što je u pretincu

Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper.