



Titration Application Pack Petrochemistry: TAN (Total Acid Number)

Proizvod #: AP0015.AT1102

EUR Cijena: Kontaktirajte nas

Nije prikazan datum isporuke

Application Pack for the determination of the Total Acid Number in Petrochem. Suited for Titrilab AT1000 Series, Automatic Titrator Models AT1102/AT1112/AT1122/AT1222.

Specifikacije

Garancija:	24 mjeseci
Kompatibilni instrumenti:	AT1102
	AT1112
	AT1122
	AT1222
Parametar:	TAN (Total Acid Number)
Primjena:	Petrochemical Applications
Što je u pretincu:	Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper, 1 L glass bottle.
Vrsta:	Application Pack

Što je u pretincu

Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper, 1 L glass bottle.