



Titration Application Pack General Chemistry: ORP I=0

Proizvod #: AP0019.AT1112
EUR Cijena: Kontaktirajte nas
Dostupno

Application Pack for ORP Endpoint Titration, I=0. Suited for Titrallab AT1000 Series, Automatic Titrator Model AT1112.

Specifikacije

Garancija:	24 mjeseci
Parametar:	ORP I=0
Što je u pretincu:	Application Pack contains Glass Syringe, Plastic Beakers, Magnetic stirring bar, Pre-programmed application on USB stick, Electrode, conical adapter, bottle stopper, 1L glass bottle

Što je u pretincu

Application Pack contains Glass Syringe, Plastic Beakers, Magnetic stirring bar, Pre-programmed application on USB stick, Electrode, conical adapter, bottle stopper, 1L glass bottle