



Titration Application Pack Generic Chemistry: ORP I>0

Proizvod #: AP0020.AT1112

EUR Cijena: Kontaktirajte nas

Otprema unutar 2 tjedna

Specifikacije

Garancija: 24 mjeseci

Parametar: ORP I>0

Što je u pretincu: Application Pack contains Glass Syringe, Plastic Beakers, Magnetic stirring bar, Pre-programmed application on USB stick, Electrode, conical adapter, bottle stopper, 1L glass bottle

Što je u pretincu

Application Pack contains Glass Syringe, Plastic Beakers, Magnetic stirring bar, Pre-programmed application on USB stick, Electrode, conical adapter, bottle stopper, 1L glass bottle