



HQ11D Digital pH meter kit, pH Gel electrode, Outdoor, 5m

Proizvod #:

HQ11D.99.105000

NEDOSTUPAN ARTIKL

Ovaj artikl više nije dostupan.

Digital Electrochemistry with HQD

Portable pH/ORP meter. Instrument with Outdoor pH Gel electrode, cable length 5m.

The digital meter/electrode system combines reliability, flexibility and ease of use. Interchangeable Intellical electrodes are recognised automatically and store all relevant data. Rugged, nearly indestructible outdoor versions with multiple cable lengths enable measurements even at previously inaccessible locations.

Replaced by NEW HQ Series.

Outstanding reliability and simple handling

Versatile electrodes for all applications, e.g. wastewater, drinking water, process water

Reliable results from inaccessible measurement locations and over long distances - even for pH

Full GLP data management

Specifikacije

AC i USB pogon:	Optional
Automatsko prepoznavanje pufera:	Yes
Funkcija zaključavanja zaslona:	Continuous / Auto-stabilization ("press to read") / At Interval
Garancija:	36 mjeseci
GLP značajke:	Date; Time; Sample ID; Operator ID
Intervali/uzbune/podsjetnik kalibracije:	Isključeno, može se odabrati od 2 sata do 7 dana
IP zaštita kućišta:	IP67
Izvoz podataka:	Download via USB connection to PC or flash memory device. Automatically transfer entire data log or as readings are taken. Power supply adapter required.
Jezici:	English, German, French, Italian, Spanish, Danish, Dutch, Polish, Portuguese, Turkish, Finnish, Czech, Russian.
Kalibracija ORP elektrode:	Predefined ORP standards (including Zobell's solution)
Komplet?:	Da
Masa:	338 g (0.75 lb) without batteries; 430 g (0.95 lb) with batteries
Mjerenje temperature:	-10 - 110 °C

Mjerna metoda:	Probe specific programmed method settings
Model:	HQ11D – pH/1 kanal
mV mjerenje :	-1500 - 1500 mV
mV mjerenje kod stabilnost očitanja:	Yes
mV rezolucija:	0.1 mV
Određene odabir pufera:	Color-coded: 4.01, 7.00, 10.01 pH; IUPAC: 1.679, 4.005, 7.000, 10.012, 12.45 DIN: 1.09, 4.65, 9.23 User-defined custom buffer sets
Parametar:	pH/Oxydo Reduction Potential (ORP)
Parameters:	pH NA NA
PC data transfer software :	HQD Series Meter Data Transfer Utility
pH Electrode calibration:	1 - 3 Calibration points Calibration summary data logged and displayed
pH Measurement:	0 - 14 pH
pH rezolucija:	Odabir između 0.001 - 0.1 pH
Pohrana podataka:	Automatic, GLP ISO compliant reading data stored with batteriescalibration details.
Pohrana podataka/memorija:	500 records/FIFO
Potvrde o sukladnosti:	CE.WEEE
Printer:	Yes, Optional
Radne poruke pogrešaka:	Clear text error messages displayed
Radni uvjeti: relativna vlažnost zraka:	90 % relative humidity (non-condensing)
Radni uvjeti: temperatura:	0 - 60 °C (32 - 140 °F)
Radno sučelje:	Soft Touch Keypad
Sastavi:	No
Sensor A:	PHC10105 (5 m rugged)
Sensor B:	NA
Sensor C:	NA
Sensors:	Sensor A: PHC10105 Sensor B: NA Sensor C:
Standardi za korisničku kalibraciju:	Yes
Što je u pretincu:	Meter package includes HQ11D Portable Meter; 4 AA batteries; PHC101 Field Rugged Gel pH Electrode, 5 m cable; Singlet pH 4 & 7 pH Buffer 10 of each; quick-start guide and user manuals.
Temperature compensation:	Automatic Temperature compensation for pH
Temperaturna rezolucija:	0,1 °C
Temperaturna točnost:	± 0.3 °C
Type of measurement:	Field Measurements
Uključuje :	Meter package includes HQ11D Portable Meter; 4 AA batteries; PHC101 Field Rugged Gel pH Electrode, 5 m cable; Singlet™ pH 4 & 7 pH Buffer 10 of each; quick-start guide and user manuals.

Ulazi:	1
Uređaj:	Portable
Vrsta sonde:	Intellical Digital SMART Standard Laboratory or Rugged Field
Vrsta zaslona:	240 x 160 pixel LCD with batteriesbacklight illumination
Zaslon:	Detailed mode/Large mode

Što je u pretincu

Meter package includes HQ11D Portable Meter; 4 AA batteries; PHC101 Field Rugged Gel pH Electrode, 5 m cable; Singlet pH 4 & 7 pH Buffer 10 of each; quick-start guide and user manuals.