



Komplet reagensa za hidrazin, 0,01-2,0 mg/L N_2H_4

Proizvod #:

LCW025

Opasno po zdravlje

EUR Cijena:

Kontaktirajte nas

Dostupno



For determination of Hydrazine by the 4-Dimethylaminobenzaldehyde method.

Specifikacije

Broj testova:	60
Čuvanje uređaja:	15 °C - 25 °C
Mjerni rang:	0,01 - 2,0 mg/L N_2H_4
Naziv metode:	4-dimetilaminobenzaldehid
Opis:	komplet reagensa za hidrazin
Parametar:	Hidrazin