



Ammonium/Potassium Combination Standard Solution: 5.000 mg/L $\text{NH}_4\text{-N}$, 20.000 mg/L K, 1000 mL

Proizvod #:

LCW892

EUR Cijena:

Kontaktirajte nas

Otprema unutar 2 tjedna

Ammonium/Potassium standard for validation purpose of ammonium process probes.

Specifikacije

Koncentracija:

5000 mg/L $\text{NH}_4\text{-N}$, 20.000 mg/L K

Parametar:

Ammonium

Uređaj:

Amtax sc