



Vrhovi za pipete 1,0-5,0 mL

Proizvod #:	LYW787
EUR Cijena:	Kontaktirajte nas
Dostupno	

Vrhovi za pipete 1,0-5,0 mL za promjenjivu pipetu 1,0-5,0 mL, fiksna pipeta od 0,2 mL i elektronička pipeta 0,2-5,0 mL

Specifikacije

Opis:	Vrhovi za pipete (1,0-5,0 mL) za BBP065, BBP164 i BBP087
Volumen/ veličina paketa:	1000 komada